Dfab

3D FOR PATIENTS

CLEAN, EASY, FAST, SMART

Irix Max & Dfab

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References:

1. Mangano FG, Yang KR, Lerner H, Porrà T, Khachatryan LG, Gordienko ID, Admakin O. 3D-printed short-span hybrid composite implant-supported restorations fabricated through tilting stereolithography: a retrospective clinical study on 85 patients with 1 year of follow-up. *J Dent.* 2024 Feb;141:104792. PMID: 38013004 DOI: 10.1016/j.jdent.2024.105095.
2. De Angelis F, D'Amario M, Jahjah A, Frascaria M, Vadini M, Sorrentino E, Biferi V, D'Arcangelo C. Flexural properties of three novel 3D-printed dental resins compared to other resin-based restorative materials. *Prosthesis.* 2024; 6(3):619-630. <https://doi.org/10.3390/prosthesis6030043>
3. Internal data, DWS R&D
4. Pontlab S.r.l. Report 1583P24 08/05/2024
5. TDS EN_Irix Max_rev.02_03/202
6. Corbani K, Hardan L, Eid R, Skienhe H, Alharbi N, Ozcan M, Salameh Z. Fracture resistance of three-unit fixed dental prostheses fabricated with milled and 3D printed composite-based materials. *J Contemp Dent Pract.* 2021 Sep 1;22(9):985-990. PMID: 35000940.
7. Mangano FG, Cianci D, Pranno N, Lerner H, Zarone F, Admakin O. True-ness, precision, time-efficiency, and cost analysis of chairside additive and subtractive versus lab-based workflows for manufacturing single crowns: An in vitro study *J Dent.* 2024 Feb;141:104792. PMID: 38013004 DOI: 10.1016/j.jdent.2023.104792.
8. IN3DLAB_G_FLD_Tecnopolimeri_REV01_V3_0922
9. Corbani K, Hardan L, Eid R, Skienhe H, Alharbi N, Ozcan M, Salameh Z. Fracture resistance of three-unit fixed dental prostheses fabricated with milled and 3D printed composite-based materials. *J Contemp Dent Pract.* 2021 Sep 1;22(9):985-990. PMID: 35000940
10. Rolando C, Rossi T, Baldi A, Comba A, Carossa M, Pera F, Alovisi M, Berrutti E, Pasqualini D, Scotti N. CAM-technique effect on cement volume and fatigue resistance of polymer-infiltrated ceramic network (PICN) restorations. Poster presented at the 2nd SIPRO (Italian Prosthodontics Society) National Congress. Florence, February 17-18, 2023 | Special Mention for Research in Prosthodontics
11. Demirel M, Donmez MB, Çakmak G, Dede DÖ, Hinz S, Yılmaz B. Effect of manufacturing trinomial and restoration thickness on the fabrication true-ness, fit, and margin quality of additively manufactured resin-based ultrathin laminate veneers *J Dent.* 2025 Feb 1:105606. DOI: 10.1016/j.jdent.2025.


Further readings

- Baldi A, Rossi T, Stura I, Comba A, Fazioni M, Rolando C, Ferrero G, Ceruti P, Scotti N. Volumetric and linear adaptation of an indirect adhesive restoration: comparison of chairside 3D printing and milling techniques. *Appl. Sci.* 2025, 15(2), 828. DOI <https://doi.org/10.3390/app15020828>.
- Corbani K, Hardan L, Eid R, Skienhe H, Ozcan M, Alharbi N, Salameh Z. Effect of material thickness on the fracture resistance and failure pattern of 3D-printed composite crowns. *Int J Comput Dent.* 2020;23(3):225-233. PMID: 32789310
- Gaetano D, Mazzone A., Baldi A., Comba A., Rossi T., Scotti N. 3-body wear of 3d printed vs milled composites: an in vitro study. DOI: 10.19256/abstract.cduo.03.2023
- Badino G., Baldi A., Fazioni M., Comba A., Rossi T., Scotti N. Anterior teeth shape modification after ortho treatment through 3d-printing technologies. DOI: 10.19256/abstract.cduo.03.2023
- IFU EN rev.03 – 01/08/2023

This document is reserved for healthcare professionals as it contains information on medical devices that may create situations of danger to the health and safety of the patient if not correctly read, understood and applied by a professional.

Dfab
3D FOR PATIENTS

For more info:
www.dfab.dental
info@rd-printing.com

J Dent. 2024 Feb;141:104792. PMID: 38013004
DOI: 10.1016/j.jdent.2024.105095.

3D-printed Short-span Hybrid Composite Implant-supported Restorations Fabricated Through Tilting Stereolithography: A Retrospective Clinical Study on 85 Patients with 1 Year of Follow-up

Francesco Guido Mangano ^a, Kyung Ran Yang ^b, Henriette Lerner ^c, Tommaso Porrà ^d, Lusine G. Khachatryan ^e, Igor Dmitrievich Gordienko ^a, Oleg Admakin ^a

a Department of Pediatric, Preventive Dentistry and Orthodontics, I. M. Sechenov First State Medical University, Moscow, Russian Federation

b MIR Dental Hospital, Daegu, South Korea

c Academic Teaching and Research Institution of Johann Wolfgang Goethe University, Frankfurt, Germany

d Private Practitioner, Timisoara, Romania

e Department of Pediatric Diseases, N. F. Filatov Clinical Institute of Children's Health, I. M. Sechenov First State Medical University, Moscow, Russian Federation

Abstract

Purpose: To report the clinical results obtained with fixed short-span (single crowns [SCs] and fixed partial prostheses [FPPs]) implant-supported hybrid composite restorations fabricated through tilting stereolithography (TSLA).

Methods: This retrospective clinical study included 85 patients who had been restored with 95 fixed short-span implant-supported hybrid composite (Irix Max[®], DWS Systems) restorations (70 SCs and 25 FPPs up to three units) fabricated with TSLA. The full-digital model-free workflow was based on intraoral implant scanning, computer-assisted design (CAD) and 3D printing using TSLA (Dfab[®], DWS Systems). The primary outcomes were the marginal adaptation, the quality of the occlusal and interproximal contact points, and the chromatic integration of the restorations, assessed independently by two experienced operators (a prosthodontist and a periodontist). A score from 1 to 5 (with 5 as the highest value, 4 for satisfactory quality, 3 for acceptable quality, and 2 and 1 as the lowest values, expressing unsatisfactory quality) was assigned by each operator to each restoration at delivery. The secondary outcomes were the survival and success of the restorations at the 1-year follow-up. The restoration was defined as successful in the absence of any complications throughout the follow-up period. A statistical analysis was conducted.

Results: For the quality of the marginal closure and occlusal and interproximal contact points, the 3D-printed hybrid composite restorations scored highly; the aesthetic integration was satisfactory. One year after placement, all restorations survived, with a low incidence (4.2% overall, 5.7% SCs) of complications (two abutment screw loosening, two decementation of the restorations, and one upper portion of the hybrid abutment decemented from the titanium base), for a success rate of 95.8%.

Conclusions: Within the limits of this study (retrospective design, follow-up limited to 1 year from the delivery, and only cemented restorations included) fixed short-span implant-supported hybrid composite crowns and bridges fabricated through TSLA were clinically precise, presenting a low incidence of complications at 1 year.

Statement of clinical relevance: The use of TSLA printing technology can open new perspectives for the treatment of small edentulous gaps with definitive implant-supported prosthetic restorations.

Copyright © 2024 Elsevier Ltd. All rights reserved.